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
English User Manual

Contents

1	Production Introduction	1
1.1	EF1800 Introduction	1
1.2	Important Safety Notes	1
1.3	Measurement Methods	2
1.3.1	Lead Wires Measurement	2
1.3.2	Measurement on Chest	5
1.3.3	Measurement on Ankle	6
1.3.4	Measurement on Hand	6
2	Unit	7
2.1	Main Unit	7
2.2	Guidance Screen	8
2.3	Utility Menu	9
2.4	ECG Waveform	9
2.5	Result Message	10
3	Operating Instructions	11
3.1	Battery	11
3.1.1	Installing Batteries	11
3.1.2	Battery Life	11
3.2	SD Memory Card	12
3.3	ECG Measurement	12
3.3.1	Notes on Taking a Measurement	12
3.3.2	Taking a Measurement	13
3.3.3	Common Causes of Inaccurate Measurement Results	17
3.3.4	Measurement Mode	17
3.4	Wave Review	17
3.5	Setting	19
3.5.1	Time/date	20
3.5.2	Language	21
3.5.3	Beeper	21
3.5.4	Default store	21
3.5.5	Format SD	21
3.5.6	Default Setting	22
3.6	Help	22
4	Trouble Shooting	22
5	Care and Maintenance	24
6	General Information about the Heart and ECG	25
6.1	General Information about the Heart	25
6.2	Electrocardiogram	26
6.3	What is Heart Rate?	28
6.4	What is Arrhythmia?	28
7	Technical Data	28
8	Product label and Mark	35
9	Customer Service	35

1 Production Introduction

1.1 EF1800 Introduction

ECG Monitor EF1800 is a portable, easy-to-use ECG monitor that can record and store the ECG measurement of your heart rhythm. Its detection principle is the same as the electrocardiograph used in hospitals, but EF1800 is portable, easy-to-use, timely detection and can auto-adjust the ECG amplitude. The EF1800 is only used for homecare. The EF1800 provides effective detecting means for heart disease patients and sub-health people to prevent heart disease. The unit can also provide the doctor with patient-related information.

- ◆ **Timely Detection:** The heart disease is more sudden and transient. When the heart is uncomfortable, people can not record the electrocardiograph as they are not in hospital. The ECG is normal while they are in hospital, so that it can not detect the cause of the disease. Thus, you need a portable ECG monitor. When you feel unwell and experience mild discomfort in the body, the EF1800 can timely detect and record the abnormality of the heart.
- ◆ **Diagnosis Assistance:** You can take a measurement by using the EF1800 any time, any where. After about 30-second measurement, the screen will display the analysis result, waveform and heart rate. The detection result is accurate and reliable, and it can assist the doctor to do pre-inspection for patient.
- ◆ **Early Prevention:** You can timely detect electrocardiogram and record the health condition of heart for early prevention, when you feel unwell.
- ◆ **Simple Operation:** 2 × AAA batteries power supply with low voltage indication, small size, light weight and portable. You can take a measurement any time, any where, and you can see the analysis result, waveform and heart rate promptly. Taking a measurement only need about 30 seconds.
- ◆ **Applicable scope:** Apply to self-monitoring and recording ECG waveform and heart rate of patient. Suitable for Heart disease patients, Sub-healthy crowd, Prevention and health care crowd (adult).

1.2 Important Safety Notes

·EF1800 is designed for monitoring the ECG waveform and knowing the heart situation. It also provides doctor with more patient-related information, in order to assist doctor to make an examination.

·EF1800 measurement results are for monitoring purposes only, self-measurement is not the same as medical diagnosis, and should never be

used as a basis for starting or modifying treatment without independent confirmation by medical examination.

·The values displayed by the unit are the values at the time of measurement. Medical conditions can change suddenly. If you notice any change in your condition, use EF1800 to take a measurement, and then consult your doctor, regardless of the measurement results.

Warning:

- Do not use this device if you have a cardiac pacemaker!
- Do not use this device with a defibrillator!
- Do not use in the presence of flammable anesthetics, drugs or pressurized oxygen (such as in a hyperbaric chamber, ultraviolet sterilizer or oxygen tent).
- Keep out of reach of small children!
- Do not expose the unit to strong shocks or vibrations, or drop or step on the unit.
- Do not apply over clothing during measurement.
- Do not take measurements when there are drops of water on your skin, such as from sweat or after bathing.
- If the unit displays “lead off” when measuring, please wipe hands or measure positions with a damp towel or disinfectant alcohol, and then repeat measurement.
- Do not take measurements where the unit will be exposing to strong electromagnetic forces or static electricity.
- Do not take measurement in a moving vehicle.
- The unit does not recommend the use of contact gel.

1.3 Measurement Methods

This unit can provide four measurement methods, the user can select one of them to take measurement according to the requirement. We recommend you to measure by the methods of lead wires measurement or measurement on chest. In the case of inconvenience, you can measure by the method of measurement on hand.

1.3.1 Lead Wires Measurement

(Lead wires are optional, and lead wires measurement must have the lead wires.)

- (1) ECG lead wires

Model: EC1000

Length: 2.1m

Color: white (+), red (—), black (N)



Function: Uses for ECG monitor EF1800

Note: Lead wires can only use PUMP production, which match the unit, model: EC1000. Else, the function and measurement results will be inaccurate.


(2) When measuring by lead wires, you can paste the disposable ECG electrodes on the right positions of the chest, you can also paste the disposable ECG electrodes on the right positions of the limbs. Put the lead wires into the lead wire socket. Put the disposable ECG electrodes on the right positions, and then connect the lead wires to the disposable ECG electrodes. (The position of the electrodes can refer to the following diagrams)

The electrode positions of the chest or the limbs:

Electrodes identifier	Color code	The electrode positions of the chest	The electrode positions of the limbs
LL (+)	White	The intersection of the left midclavicular line and the fifth rib	Inner of the left leg
RA (-)	Red	The intersection of the right midclavicular line and the second rib	Inner of the right arm
RL (N)	Black	The intersection of the right midclavicular line and the fifth rib	Inner of the right leg

Positions of the chest



positions of the limbs

Note:

- The hands can not touch the metal electrodes on the both sides when measuring by lead wires, otherwise it will interfere the measurement.
- If the unit displays “lead off” when measuring, please wipe the measure positions with a damp towel or disinfectant alcohol, and check the lead wires, and then repeat measurement.

(3) Use method of the disposable ECG electrode:

Lead wires should be used with the ECG electrode which is disposable, when measuring by lead wires. Disposable ECG electrode has the advantages of easy to use, paste firmly, not easy to happen the cross-infection and so on. But in order to ensure that the unit obtains better ECG information when use it, attention should be paid to the following problems.

- Should keep the paste position cleaning before using the disposable ECG electrode, user should wipe the paste position with the cotton swab or disinfectant alcohol. If the skin is dirty or the stratum corneum is too thick, it may lead to resistance increases, current decreases, and then the occurrence of adverse conduction.

- Accurately paste the electrodes on the part, paste should be gently pressed, to ensure that the conductive pressure-sensitive adhesive in the center of the electrodes close to the skin.

- After the disposable ECG electrodes are torn out from the skin, it can not be used again. Otherwise, it may cause the measurement result is inaccurate or the ECG information can not be detected.


- The time of disposable ECG electrode pasting on the skin should not exceed 24 hours. If the time exceeds the period, user should tear out or replace the disposable ECG electrode in time, so as to avoid the detected signal is adverse or the skin allergy.

Note: When you want to buy the disposable ECG electrode, please go to the normal pharmacy (it must also sell the medical equipment) or manufactures. Please pay your attention to its manufacturing license and product registration when buying the disposable ECG electrode.


1.3.2 Measurement on Chest

Hold the unit with your right hand. Put right index finger on the two right electrodes, and place the left electrode on the intersection of the left midclavicular line and the 5th rib. The position of the left electrode can refer to the follow diagrams. If you are unsure, please consult your doctor.

Male: Place the left electrode on bare skin about 5 cm below your left nipple.

Female: Locate the EF1800 at the lower end of the breastbone and move horizontally to the left to the middle of the part of the chest, if necessary, gently lift the breast and place the left electrode under the left breast. If the left electrode comes in contact with the brassiere, please remove it, and make sure the electrode contact with the breast directly. If you are unsure, please consult your doctor.

The intersection of the left midclavicular line and the 5th rib.

Note: If the unit displays “lead off” when measuring, please wipe hand or measure positions with a damp towel or disinfectant alcohol, and then repeat measurement.

1.3.3 Measurement on Ankle

Hold the unit with your right hand. Put right index finger on the two right electrodes, and place the left electrode on the inner of the left ankle.

Note: If the unit displays “lead off” when measuring, please wipe hand or measure positions with a damp towel or disinfectant alcohol, and then repeat measurement.

1.3.4 Measurement on Hand

Hold the unit with your right hand. Put right index finger on the two right electrodes, and place the left electrode on your left palm.


Note:

- The hands should not touch each other when measuring.
- If the unit displays “lead off” when measuring, please wipe hands with a damp towel or disinfectant alcohol, and then repeat measurement.

2 Unit


2.1 Main Unit


- ①. **Right electrodes:** Be composed of two metal electrodes. Place your right index finger on these when taking a measurement.
- ②. **Left electrode:** Place this on your left hand or on the inner of the left ankle, or against your chest when taking a measurement.
- ③. **TFT display:** Displays the wave pattern after measurement, and displays menus and setting options.
- ④. **Battery cover:** Open this cover to insert batteries.
- ⑤. **Lead wire socket:** Put the lead wires into the socket when using lead wires measurement method.
- ⑥. **(ON/OFF) button:** Turn on/off the power.
- ⑦. **START button:** Starts measurement.
- ⑧. **MENU button:** Displays or closes the menu screen.
- ⑨. **↶↷ button:** Moves the desired menu or wave pattern to the left or up and increases the number of the selected field.
- ⑩. **↵↘ button:** Moves the desired menu or wave pattern to the right or down and decreases the number of the selected field.
- ⑪. **OK button:** Selects the selected menu option or ECG data.
- ⑫. **SD card cover:** open this cover to insert the SD memory card.

2.2 Guidance Screen

When turn on the power button, the follow screen is displayed. It displays a standard ECG wave, company logo and company name.

Then, the screen displays a reminder that self-measurement is not the same as medical diagnosis.


After the message, the guidance screen is displayed, it used to guide user how to take a measure.

2.3 Utility Menu

TFT screen is used to display the guidance screens and results of ECG readings, adjust settings, and manage ECG data.

Turn on the power, and then press the MENU button back to the utility menu. This unit contains 4 options:

ECG measurement: While maintaining the right measurement posture, press the START button. During measurement, the screen displays the ECG waveform and the progress bar at the bottom of the screen displays measurement progress. About 30 seconds later, the heart rate and analysis result are displayed. (Refer to 3.3 ECG measurement)


Wave Review: This allows you to view ECG waves and results stored in the memory. (Refer to 3.4 Wave Review)

Setting: The setting menu allows you to set the time/date, and adjust language, beeper or default setting. (Refer to 3.5 setting)

Help: This displays the guidance screens that show you how to take a measurement. (Refer to 3.6 help)

2.4 ECG Waveform


The screen shows as the following diagram when taking a measurement. It contains ECG waveform, time and the progress bar at the bottom of the screen displays measurement progress.

- ①. ECG waveform.
- ②. The progress bar.
- ③. Calibration guide: Serves as a guide when your doctor analyses the stored ECG data.
- ④. Battery level indicator.
- ⑤. Magnification: contains: 1/2, 1, 2, 3.
- ⑥. Indicates a measurement is taking.
- ⑦. Time.

2.5 Result Message

A measurement result message is displayed when you review the ECG data. It displays the heart rate, the result and the relevant operation.

3 Operating Instructions

3.1 Battery

3.1.1 Installing Batteries

- (1) Slide the battery cover in the direction of the arrow on battery cover, and lift the end of the cover upwards.
- (2) Insert two LR03(AAA) batteries as indicated. Make sure that the polarity of the batteries is correct.

- (3) Close the battery cover, then while gently pressing the cover down, slide it in the direction against the arrow on battery cover.

3.1.2 Battery Life

- (1) When the “low battery!” message appears, replace both batteries with new ones.
 - Always turn off the unit before replacing the batteries.
 - Dispose of the used batteries according to applicable local regulations.
- (2) Battery life

The battery life can be confirmed in the upper right of the display.

The states of battery life:

—— Full power.

—— Half remaining.

—— Remaining power is low. Replace the batteries with new ones.

Notes:

- Do not use new and used batteries together.
- Do not use batteries of different types together.
- Do not insert the batteries with the polarities in the wrong direction.
- Do not dispose of batteries in fire.
- Remove the batteries from this unit when you are not going to use it for a long period of time (approximately one month or more).
- Please use the batteries which produced by the regular manufactures, and they must conform to national standards.

3.2 SD Memory Card

EF1800 includes an SD memory card for storing ECG data, which can then be viewed by your doctor later.

(1) insert the SD memory card

Check that the unit is turned off, open the SD memory card slot cover and insert the SD memory card, make sure that the card is correctly aligned, and then press the SD memory card down gently, the SD memory card is fixed in the slot, finally close the SD memory card slot cover.

(2) eject the SD memory card

Check that the unit is turned off, open the SD memory card slot cover, press the SD memory card gently; the lock is released and the card is ejected slightly, and then remove the SD memory card, finally close the SD memory card slot cover.

Notes:

- Never insert or eject an SD memory card while the unit is turned on.
- Do not store any data on the card from other devices other than this unit.
- Make sure that the write-protect switch is not locked, before inserting the SD memory card. If the write-protect switch is locked, you will not be able to store new ECG data on the SD memory card.
- When inserting the SD memory card, do not try to force the card into the slot. If the SD memory card can not be inserted into the slot, check whether the direction is correct, otherwise this could damage the card or the unit.

3.3 ECG Measurement

3.3.1 Notes on Taking a Measurement

In order to obtain accurate ECG information, it is important that measurements are taken correctly. Before taking a measurement for the first time, read these instructions carefully and make sure you observe the following points each time you take a measurement.

- Make sure that the electrodes are directly touching your skin.
- If your hands or skin are dry, wipe them with a damp towel so that they are slightly moist.
- If the unit displays “lead off” when measuring, please wipe hands or measure positions with a damp towel or disinfectant alcohol, and then repeat measurement.
- If the electrodes are dirty, wipe any dirt off with a soft cloth moistened with disinfectant alcohol or a cotton swab.
 - Take care not to use too much disinfectant alcohol when cleaning the electrodes.
 - Do not wipe or get disinfectant alcohol on parts other than the electrodes.
- Keep quiet during measurement. Movement, including talking, coughing, or sneezing can affect the measurement results during measurement.

·If the position of your finger or arm is inappropriate during measurement, the results may be unreliable, so be sure to check their position before taking a measurement. (Refer to 1.3).

·If the signal is not ideal, please use other measurement methods.

3.3.2 Taking a Measurement

3.3.2.1 Measuring Use Electrodes

- (1) Press the  button more than 3 seconds to turn on the unit.
The unit displays a reminder that self-measurement is not the same as medical diagnosis and the guidance screen.
- (2) Make sure that you are in a relaxed position and seated comfortably with your back straight. You can also take a measurement while you are lying.
- (3) Fit your right index finger closely on the right electrodes.
- (4) Place the left electrode on the left hand or the inner of the left ankle or bare skin below your left nipple. (refer to 1.3) .
- (5) Push START button and hold stable for about 30 seconds.

Note: Take care not to press the START button until you are ready to take the measurement.



3.3.2.2 Measuring by Lead Wires

(Lead wires are optional, and lead wires measurement must have the lead wires.)


(1) Put the lead wires into the lead wire socket. Put the disposable ECG electrodes on the right positions, and then connect the lead wires to the disposable ECG electrodes.

The electrode positions of the chest or the limbs:

Electrodes identifier	Color code	The electrode positions of the chest	The electrode positions of the limbs
LL (+)	White	The intersection of the left midclavicular line and the fifth rib	Inner of the left leg
RA (-)	Red	The intersection of the right midclavicular line and the second rib	Inner of the right arm
RL (N)	Black	The intersection of the right midclavicular line and the fifth rib	Inner of the right leg

Positions of the chest

positions of the limbs

Notes:

- The position of the disposable ECG electrodes can refer to the above diagrams.
- Lead wires measurement must have the lead wires, and lead wires can only use PUMP production which matches the unit, model: EC1000.
- Lead wires should be used with the ECG electrode which is disposable, use method of the disposable ECG electrode can refer to 1.3.1(3).

- (2) Press the  button more than 3 seconds to turn on the unit.


The unit displays a reminder that self-measurement is not the same as medical diagnosis and the guidance screen.

- (3) Make sure that you are in a relaxed position, and you are seated comfortably with your back straight. You can also take a measurement while you are lying.
- (4) Push START button and hold stable for about 30 seconds.


The measurement takes about 30 seconds to complete. After press the START button, there are 3 seconds for preparation, after 3 seconds it start to measure. During measurement, the screen displays the ECG waveform and the progress bar at the bottom of the screen displays measurement progress.

If the contact of the electrodes becomes loose, the screen displays “Lead off” at the middle of the top.

A measurement result message is displayed when the measurement has been completed. The message includes result, heart rate, and asking return or save.

If you press the MENU button while the result message is displayed, it will return to previous menu, and the ECG data will not be save. If you press the OK button while the result message is displayed, the ECG data will be saved in default memory automatically. Default store setting refers to 3.5.4. If the unit's memory is set to the default memory, and the unit's memory is full, the ECG


data will be saved in SD memory card automatically.

Notes: ·If there is no SD card in the unit, the ECG data will be saved in unit's memory automatically.

·The unit memory can store 40 measurements, and the SD memory card which we offer can store 1G data.

If you not press the OK button or MENU button while the result message is displayed. The unit will save the ECG data automatically after about 15 seconds.

If the SD memory card and unit's memory are full, or there is no SD card in the unit, and the unit's memory is full, the screen display a reminder that memory full! Overwrite record?

When the remind message is displayed, you can press  button or  button to select "Yes" or "No", and then press OK button to confirm.

When you select "Yes", the latest measurement data will overwrite the first measurement data. When you select "No", the latest measurement data will not be stored.

Notes:

·When the screen displays "Lead off" at the middle of the top, please check the positions of the electrodes. If the positions of the electrodes are correct, it may be caused by the skin which is too dry, at this time please wipe hands or measure positions with a damp towel or disinfectant alcohol, and then repeat measurement.

·If the "low battery!" message is displayed immediately after measurement, and the unit is shut off automatically before the data is stored, the ECG data may not be stored. If this happens, please repeat the measurement after replacing the batteries.

·If you forget to turn the unit off, it will automatically shut itself off after 74 seconds.

3.3.3 Common Causes of Inaccurate Measurement Results

- (1) The index finger is not placed closely on the electrode.
- (2) The unit is held incorrectly.
- (3) Unit placed on clothing during measurement.
- (4) The electrode is placed incorrectly:
 - a. Put the left finger on the right electrodes.
 - b. The left electrode is against the finger and the right electrodes are against the left palm, left ankle or chest.
- (5) The electrodes of the lead wires are placed incorrectly.

3.3.4 Measurement Mode

There are 9 measurement modes regarding the measurement analysis.

1	Normal
2	Beat Fast
3	Beat Slowly
4	Suspected premature beat
5	Suspected ve_bigeminy
6	Suspected ve_trigeminy
7	Suspected missed_beat
8	Irregular heart rate
9	Signal bad measure again

Important:

· Do not interpret measurement mode yourself, self-measurement is not the same as medical diagnosis, please consult your doctor.

3.4 Wave Review

The Wave Review is used to display the list of ECG data stored in the unit memory or in the SD memory card.

(1) Press the MENU button while the unit is turned on.

(2) Select the wave review icon, the Wave Review is selected.



(3) Then press “OK” button, you can choose “device” or “SD card” to review the stored data.

Note: if there is no SD card in the unit, you can not choose to review the data in the SD card.

(4) If you choose “device” or “SD card” to review the stored data, and then press the OK button, the list of ECG data that stored in device or SD card is displayed. If you review the data in the “device”,


the list shows as the right figure. If you review the data in the “SD card”, there has a SD symbol in the list. Press \uparrow/\leftarrow button or \downarrow/\rightarrow button to select the ECG data you wish to view. Each stored ECG data is displayed by date and time.


(5) Select the ECG data you wish to view from the records, then press the OK button to display waveform. Press \uparrow/\leftarrow button or \downarrow/\rightarrow button to view the detail of the whole 30 second measurement.

(6) Press the MENU button again to view the measurement comments and heart rate. The results review in the “device” and “SD card” are not same, they are showed in the following:

Device

SD card


(7) If you choose to review the result in device, you can operate: Copy all, delete, delete all and return. If you choose to review the result in SD card, you only can operate: delete, delete all and return, while you can not copy data.

Select “Return”, then press OK button to return to ECG lists.

If you select “Delete”, then press OK button, it reminds that delete this ECG data. Select “Yes”, and press OK button, the ECG data will be deleted.

If you select “Delete All”, then press OK button, all the ECG data stored in the unit or SD card will be deleted.

If you select “Copy All”, then press OK button, the reminding will be displayed, it showed as the right figure. Choose “Yes”, and then press OK button, the all data in the device will be copied to SD card.

Notes:


- The most recent measurement result is displayed at the top of the list and is automatically highlighted, with any previous measurement results listed below it in chronological order.

Press button to move up, or press button to move down.

- When you delete all data in the SD card in EF1800, you may take longer, please be patient. If you want to delete all data in SD card promptly, you must delete them in the computer.
- There will remind that “copying No. n record” during copy, at that time, please not to remove SD card, if you remove SD card during copy, the data may be damaged. If the unit displays the messages “low batteries, operation abort” during copying, you need to change the batteries, and then copy them again.
- If there is no SD card in the unit during copying, there will display message “SD is write-protect”, you can press MENU button to return.

3.5 Setting

You can use the setting to set the time/date, and adjust various settings.

Basic Operation

The following procedure describes how to select and adjust a setting. Refer to the following pages for details of the options available for each setting.

- (1) Press the MENU button, while the unit is on, to display the menu screen..
- (2) Press the  button or  button to select Setting.
- (3) Press the OK button. Initially, the Time/date setting is selected.
- (4) Press the  button or  button to select the desired setting, then press OK button to see the options for that setting.
- (5) Adjust the settings then press the OK button to confirm the changes.

3.5.1 Time/date

Always set the Time/date before using the unit for the first time. Check the Time/date is correct before using the unit, reset them if necessary. The Time/date is important indicators of when a measurement was taken.

- (1) Operating according to the Basic Operation until the Time/date is displayed.
- (2) Press the  button or  button to adjust “day”.

Press the  button to increase the number.

Press the  button to decrease the number.


- (3) Press the OK button to confirm the day setting.
The day setting is set and the month setting is highlighted.
The display cycles through the settings:
Day>month>year>hour>minute>second

- (4) In the same way, press the  button or  button set the month, year, hour , minute, and second

After you complete the Time/date setting, press OK button to confirm. Press MENU button to return to the menu screen at any time.


3.5.2 Language

This setting specifies the language to be used for the screen display. Press the \uparrow/\leftarrow button or \downarrow/\rightarrow button to select the desired language, then press the OK button to confirm your selection.

3.5.3 Beeper

Press the \uparrow/\leftarrow button or \downarrow/\rightarrow button to turn the beeper on or off, then press the OK button to confirm your selection.

3.5.4 Default store

After entering the default store menu, press the \uparrow/\leftarrow button or \downarrow/\rightarrow button to select “device” or “SD card”, and then press the OK button to confirm your selection. If you select “device”, the data will be saved in the unit’s memory automatically after measurement. If you select “SD card”, the data will be saved in the SD card automatically after measurement.


3.5.5 Format SD


After entering the format SD menu, asked whether “sure to delete all?” Press

the  button or  button to select “Yes”, and then press the OK button to confirm your selection, all data in the SD card is deleted. If you select “No”, the data will not be formatted.

Note: format SD card will not affect the data in the device.

3.5.6 Default Setting

Press the  button or  button to select “Device” or “SD Card”, and then press the OK button to confirm your selection.


Notes:

- If you press the OK button while “Yes” is selected, all of the ECG data in the device will be formatted, while the data in the SD card are not affected, the language will come back Chinese and the default memory will come back to device.
- If you press the OK button while “No” is selected, the unit will not have any changes.

3.6 Help

This help menu can guide user how to take a measure.

Press the MENU button, while the unit is on, to display the menu screen.

- (1) Press the  button or  button to select Help.
- (2) Press the OK button, then the guidance screen is displayed.
- (3) Press the  button or  button to view more help information.
- (4) Press MENU button to return to the menu screen.

4 Trouble Shooting

User may encounter some problems during using the EF1800, this moment user can solve the problem referring to the following table.

Problem	Cause	Solution
Nothing happens when you press the  button more than 3 seconds.	Not have batteries.	Insert batteries with new ones correctly.
	The batteries are worn out.	Replace both batteries with new ones.
	The batteries are inserted incorrectly.	Insert the batteries with their polarities correctly aligned.
Measurement finishes but the message “Signal bad, measure again” appears.	An ECG wave-pattern could not be detected.	Check the user manual and repeat measurement.
	Electrode contact was broken before 30 seconds of data has been saved.	Repeat measurement, and make sure that the measurement parts maintain contact with the electrodes for the full 30 seconds of the measurement.
Measurement starts but the message “Lead off” appears.	The electrodes are not making good contact with your bare skin.	Check the user manual and repeat measurement.
	Your body is too tense.	Relax and repeat the measurement.
	Your skin is dry.	Moisten your skin with a damp towel.
	The electrodes have been applied over clothing.	Apply the electrodes to bare skin.
	The ECG signal is too weak.	Make sure that the electrodes are applied correctly. However, for some individuals measurement may not be possible.
	Lead wires not fully insert into the lead socket when measuring by lead wires.	Check the lead wires, and then repeat measurement.
The measurement signal confusion	Every body's skin is different, as the stratum corneum is too thick, leads to resistance increases and so on, and may affect the measurement.	If the signal confusion, when user measurement on hand, suggest to measure on left ankle; if the signal still confusion when measurement on ankle, suggest to measure on chest.

Problem	Cause	Solution
	The unit is exposing to strong electromagnetic forces.	Check whether there is electromagnetic interference, for example, somebody is using cell phone around user, etc. If there have the electromagnetic interference, turn off the equipment which may occur electromagnetic interference, or repeat measurement after change the environment.
The measurement signal is too weak when measuring by lead wires.	The disposable ECG electrodes are used again.	The disposable ECG electrode can not be reuse. Replace new disposable ECG electrodes, and then repeat measurement.
	The positions of the electrodes are inaccuracy.	Check the user manual and repeat measurement.
The SD memory card in the unit, but the unit can not recognize it.	The SD memory card can not be recognized.	Replace a new SD memory card.
The SD memory card can not be used.	The SD memory card may be damaged.	Format the SD memory card, or use a new SD memory card.
ECG data has disappeared	Device error.	Please contact PUMP customer service.
Other problems		Turn off the unit, then turn on the unit and try again. If the date and time or other settings have been reset to their default values, set them again.

If the problem still occurs, or happen to some other problems, please contact PUMP customer service.

5 Care and Maintenance

The EF1800 should be maintenance at a period, in order to prolong life of the EF1800 and its each component. Specific requirements are as follows:

- (1) Keep the electrodes cleaning after use the unit. If the electrodes are dirty, wipe

any dirt off with a soft cloth moistened with disinfectant alcohol or a cotton swab.

- (2) Do not wipe or get disinfectant alcohol on parts other than the electrodes.
- (3) Do not expose the unit to high temperatures, direct sunlight, high humidity, dusty locations or strong vibration.
- (4) Do not disassemble, repair, or modify the unit. If you have any question, please connect the customer service department of the manufacturer.

User should do the regular safety check for ECG monitor except the maintenance, to prevent accidents. Regular safety check should be doing as the following:


- (1) Label ---- check whether the label is clear.
- (2) Lead wires ---- check whether the wire surface of the lead wires have damage, scratch, bamboo and other phenomena. Whether the clips of lead wires have fracture etc.
- (3) SD card ---- check whether the SD card is convenient to insert or eject, whether the function of the SD card is normal.

6 General Information about the Heart and ECG

6.1 General Information about the Heart

The heart is an important organ of the human circulatory system. Relying on the rhythmic contraction and relaxation of the heart, blood can be able to flow incessantly in a closed circulation system, and oxygen transported to various parts of the body tissues, carbon dioxide emissions from the human body, so that life can be maintained

In order to analyze the principle of ECG, first of all we need to explain the process of heart activity. Electrical conduction system of the heart shown below:


The normal electrical conduction in the heart allows the impulse that is generated by the sinoatrial node (SA node) of the heart to be propagated to the atria and

ventricles, cause the excitement of the whole heart. Specifically, the SA node generates an impulse, and be propagated to the right and left atria, causing the right and left atria to contract and pump blood into the right and left ventricles respectively. As the impulse from the atrium to the ventricle has a special pathway, the impulse arrived to the ventricle in a short time, causing the ventricle to contract. Therefore, in each cardiac cycle, the various parts of the heart occur the electrical signal during the excitement, the changes of the signal in direction, way, sequence and timing has certain regularity. The changes of the electrical signal spread across the surrounding body electrolyte of the heart to body surface, so that the electrical signal various parts of the body in each cardiac cycle also have regular changes.

6.2 Electrocardiogram

An electrocardiogram (ECG) is the representation of the electrical activity of the heart as it is recorded from the body surface. It indicates the electrical activity of the heart during the excitement generated, conduction and the restoration. Typical ECG waveform shown below:

(1) Typical ECG Wave:

P wave: The P Wave indicates the left and right atria depolarization. The normal electrical conduction in the heart allows the impulse that is generated by the sinoatrial node of the heart to be propagated to the atria first. So the first deflection of the waveform is a small upward wave called the P wave. The initial portion of the P wave is largely a reflection of right atrium depolarization and the terminal portion is largely a reflection of left atrium depolarization. P wave is small and obtuse; it has slight difference in each lead. The P wave should be no larger than 0.11s, and the

maximum amplitude no larger than 2.5mm.

QRS complex: The QRS complex indicates the left and right ventricular depolarization. Typical QRS complex contain three closely link wave, initial downward deflection after the P wave called the Q wave, the first upward deflection after the P wave called the R wave, and the downward deflection after R wave called the S wave. The width of the QRS complex called QRS time limitation, it indicates the required time of ventricular depolarization. The normal QRS complex should be no larger than 0.10s.

T wave: The T wave is representation of ventricular repolarization. The amplitude of the normal T wave is about 0.1~0.8mV, Higher in the R wave, T wave amplitude should not be less than 1/10 of R wave.

U wave: The U wave located after the T wave. The U wave may be raise with hypokalemia and ventricular hypertrophy. The U wave inversion may represent hyperkalemia.

(2) Typical Interval and Typical Segment

P-R interval: The PR interval is measured from the beginning of the P wave to the beginning of the QRS complex. It represents the time required for an electrical impulse to depolarize the atria and reach the electrical conduction system of lower part of the heart, the ventricular. This period has the trend of lengthening with the growth with age.

QRS interval: The QRS interval is measured from the beginning of the R(Q) wave to the end of the S wave. It represents ventricular depolarization.

Q-T interval: The Q-T interval is measured from the beginning of the Q wave to the end of the T wave. It represents the total time of ventricular depolarization and repolarization, usually less than 0.4s, it is easily influenced by the heart rate.

PR segment: The PR segment is measured from the terminal portion of the P wave to the beginning of the QRS complex. The normal PR segment is close to the baseline.

ST segment: The ST segment is measured from the terminal of the QRS complex to the beginning of the T waves. It represents the repolarization of the ventricle. The normal ST segment is also close to the baseline, and the distance between the baseline and the ST segment is less than 0.05mm.

6.3 What is Heart Rate?

Heart rate (commonly called the number of the heartbeat) means the heart throb in the unit of time, commonly means the number of the heartbeat per minute. The normal adult heart beats at about 60 BPM (beats per minute) to 100BPM, when resting. The average heart rate is about 75 BPM. In adults, the heart rate in females is faster than in males. One's heart rate will slow down when resting or sleeping, and when movements or emotionality, the heart rate will accelerate. In certain drugs or neurohumoral factors, the heart rate will speed up or slow down. People who is regular physical exercise, the normal heart rate also slowed down, which is a natural physiological phenomenon.

6.4 What is Arrhythmia?

Arrhythmia is a condition where the heartbeat rhythm is abnormal due to flaws in the bio-electrical system that drives the heartbeat. Arrhythmia can occur in patients with heart disease, also can occur in normal people.

The normal rhythm is stabile and regular in accordance with some certain rhythm, the normal adult heart beats at about 60 BPM to 100BPM. In regular situation, there generated an impulse first before the heart beat. People can not feel uncomfortable about the regular heartbeat in general, but if the heart happens to arrhythmia, people will feel unwell, such as feeling tame, nervous, dizziness, and fatigue and so on. Different people have different feelings, because arrhythmia also has many different types.

7 Technical Data

EF1800 technical data:

Heart Rate Range: 30~220BPM

Heart Rate Measurement Errors: $\leq \pm 3\%$.

Display: TFT color display

Memory: the unit memory can store 40 measurements
the SD memory card can store 1G data.

Electrical Safety: Internally powered equipment, BF-Type device

Power Supply: 2 LR03(AAA) batteries

Operating Temperature and Humidity: Temperature: $+5^{\circ}\text{C} \sim +40^{\circ}\text{C}$, Humidity: $<80\%$

Storage Temperature and Humidity: Temperature: $-20^{\circ}\text{C} \sim +60^{\circ}\text{C}$, Humidity: $10\% \sim 95\%$

Range of operating atmospheric pressure: 70.0kPa ~ 106.0kPa

Range of storage

atmospheric pressure:	50.0kPa ~ 107.4kPa
Weight:	Approximately 100g
Dimensions:	125mm×68mm×20mm
Attachments:	Batteries 2×LR03(AAA), SD memory card (1G), User manual, Registration card, Guarantee card, Customer service card
Options:	ECG lead wires EC1000

= Type BF

**Notes:**

- Specifications may be changed without prior notice.
- Disposal of this product and used batteries should be carried out in accordance with the national regulations for the disposal of electronic products.

Important information regarding Electro Magnetic Compatibility (EMC)

With the increased number of electronic devices such as PC's and mobile (cellular) telephones, medical devices in use may be susceptible to electromagnetic interference from other devices. Electromagnetic interference may result in incorrect operation of the medical device and create a potentially unsafe situation.

In order to regulate the requirements for EMC (Electro Magnetic Compatibility) with the aim to prevent unsafe product situation, the IEC60601-1-2 standard has been implemented. This standard defines the levels of immunity to electromagnetic interferences as well as maximum levels of electromagnetic emissions for medical devices.

This medical device manufactured by pump conforms to this IEC60601-1-2:2007 standard for both immunity and emissions.

Nevertheless, special precautions need to be observed:

- Do not use mobile (cellular) telephones and other devices, which generate strong electrical or electromagnetic fields, near the medical device. This may result in incorrect operation of the unit and create a potentially unsafe situation.

Recommendation is to keep a minimum distance of 7 m. Verify correct operation of the device in case the distance is shorter.

Further documentation in accordance with IEC60601-1-2:2007 is available at pump at the address mentioned in this user manual.

Guidance and manufacturer's declaration**Guidance and manufacturer's declaration – electromagnetic emissions**

The [EQUIPMENT or SYSTEM] is intended for use in the electromagnetic environment specified below. The customer or the user of the [EQUIPMENT or SYSTEM] should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
-----------------------	-------------------	---

RF emissions CISPR 11	Group 1	The [EQUIPMENT or SYSTEM] uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable	

Guidance and manufacturer's declaration – electromagnetic immunity

The [EQUIPMENT or SYSTEM] is intended for use in the electromagnetic environment specified below. The customer or the user of the [EQUIPMENT or SYSTEM] should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. If ESD interfere with the operation of equipment, counter measurements such as wrist strap, grounding shall be considered.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.

Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles < 5% UT (>95% dip in UT) for 5 sec	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the [equipment or system] requires continued operation during power mains interruptions, it is recommended that the [equipment or system] be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

**Guidance and manufacturer's declaration – electromagnetic immunity –
for EQUIPMENT and SYSTEMS that are not LIFE – SUPPORTING**

Guidance and manufacturer's declaration – electromagnetic immunity			
The [EQUIPMENT or SYSTEM] is intended for use in the electromagnetic environment specified below. The customer or the user of the [EQUIPMENT or SYSTEM] should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3V 3V/m	Portable and mobile RF communications equipment should be used no closer to any part of the [EQUIPMENT or SYSTEM], including cables, than the recommended
Radiated RF IEC	3 V/m 80 MHz to 2.5		

61000-4-3	GHz		<p>separation distance calculated from the equation applicable to the frequency or the transmitter.</p> <p>Recommended separation distance</p> $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2.3 \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
-----------	-----	--	--

Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM –
For EQUIPMENT and SYSTEMS that are not LIFE – SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and the [EQUIPMENT or SYSTEM]

The [EQUIPMENT or SYSTEM] is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the [EQUIPMENT or SYSTEM] can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the [EQUIPMENT or SYSTEM] as recommended below, according to the maximum output power of the communications equipment.


Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.16 \sqrt{P}$	80 MHz to 800 MHz $d = 1.16 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.33 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

8 Product label and Mark

Notes:

Symbol for "date of manufacture"

Symbol for "manufacturer"

9 Customer Service

Product Name: ECG monitor

Model: EF1800

Manufacturer: **SHENZHEN PUMP MEDICAL SYSTEM CO., LTD.**


Address: 2/F West, M-7 Sinosteel Building, Maqueling Estate, Hi-Tech Industrial Park, Nanshan District, Shenzhen 518057, China

Tel: 86 755 26710795, 26710793

Service line: 86 4006 755 009 Fax: 86 755 26012025

E-MAIL: service@bpump.com.cn

Website: <http://www.bpump.com.cn>

Shanghai International Trading Corp. GmbH (Hamburg)

Address: Eiffestrasse 80, 20537 Hamburg, Germany